

SEROLOGY MULTIMARKER RANGE



QConnect ToRCH G



Storage and Stability

Sealed: 24 months at 2 to 8°C
 After opening: 90 days at 2 to 8°C

Analytes included

Anti-Rubella IgG
 Anti-Toxoplasma IgG
 Anti-CMV IgG
 Anti-EBV VCA IgG
 Anti-EBV EBNA IgG
 Anti-HSV 1 IgG
 Anti-HSV 2 IgG

CAT. NO.	CONTENT	KIT	
SR03013	Positive control	4 x 2.5 mL	4 x 2.5 mL
SR03014	Negative control Positive control	2 x 2.5 mL 2 x 2.5 mL	4 x 2.5 mL

For further information see instructions for use.



QConnect ToRCH M



Storage and Stability

Sealed: 24 months at 2 to 8°C
 After opening: 90 days at 2 to 8°C

Analytes included

Anti-Rubella IgM
 Anti-Toxoplasma IgM
 Anti-CMV IgM
 Anti-EBV VCA IgM
 Anti-HSV 1 IgM
 Anti-HSV 2 IgM

CAT. NO.	CONTENT	KIT	
SR03023	Positive control	4 x 2.5 mL	4 x 2.5 mL
SR03024	Negative control Positive control	2 x 2.5 mL 2 x 2.5 mL	4 x 2.5 mL

For further information see instructions for use.